

In This Issue

Product Information Helfpul Hints / FAQs Accessories

# **MVE Chart Tech Tips**

#### **PRODUCT INFORMATION**

#### Updated TEC Connect Program (Ver. 1.1.0.0)

An updated version of the TEC Connect program is available here.

Additions/Changes:

- Added Settings Manager functionality for Vario Pro controller
- Improved functionality with the TEC 2000
- Added Excel 2007, 2010 and 2013 compatibility to the download function. Graph will be automatically generated for data when downloaded.
- Added Fill Start, Fill Stop to be compatible with TEC 2000 and TEC 3000
- Added Alarm Mute button to be compatible with all models
- Added Alarm Monitoring for the TEC 2000/3000
- Added Driver Installation Feature to the File menu

To update to the new version, follow the steps below:

- Before downloading version 1.1.0.0 of TEC Connect, make sure to uninstall any previous versions of TEC Connect. (If using Windows 7, go to Control Panel > Programs > Programs and Features. Right click the TEC Connect program and select Uninstall.)
- 2. Once any old versions of TEC Connect have been uninstalled, download the TEC Connect Software package from <u>here.</u>





March 2015

- 3. Extract all files from the downloaded folder.
- 4. Run the TEC\_Connect\_Setup application.
- 5. TEC Connect should open after installation. If not, navigate to the program in the file browser and open it.

Note: The TEC COM cable drivers are downloaded along with the TEC Connect program, but if the cable drivers have been installed previously, there is no need to reinstall the drivers. If this is the first time TEC Connect is being installed on a computer, the TEC COM cable drivers will need to be installed also. The TEC Connect program can check if the cable drivers are installed by selecting **File > Cable Drivers Check/Install.** The cable drivers can also be installed manually using the **\_install.bat** file in the TEC COM Cable Drivers folder.

File Help		
Settings	mation	
Cable Dri	vers Check/Install	
Close		
Connect	Response	
Unit ID: V Auto Fin Single Control	Controls TEMP A TEMP B Start Fill Stop Fill	LN2 LM Change Settings
Download Data Range	36000 A	Browse

#### CryoCart Replacement Casters

Chart offers replacement casters for the CryoCart. The design of the CryoCarts have changed over the years, so it is important to obtain the part number and serial number of the CryoCart in order to select the proper repalcement casters, as the casters between the two models are not interchangeable. Two of each type of caster listed will outfit one CryoCart. Part numbers 14448703 and W3170LT are the current model CryoCarts. The replacement casters for these models are:

PN 46-912042 - 4" non-pneumatic, swivel caster with brake PN 46-912043 - 4" non-pneumatic, non-swivel caster

Part number 12966911 is the previous model CryoCart. The standard casters on this model were pneumatic:

PN 12967219 - 6" pneumatic, swivel caster with brake PN 12967227 - 6" pneumatic, non-swivel caster

Alternatively, there are non-pneumatic casters available for this model, listed below:

PN 3110599 - 5" non-pneumatic, swivel caster with brake PN 3110609 - 5" non-pneumatic, non-swivel caster

### HELPFUL HINTS / FAQ's

Q: How can I replace the lid liner on an XLC 1211 or MVE 1223 model freezer?

A: The Lid liner originally used on the XLC1211 and MVE1223 is obsolete. The liner with a cork insert was replaced by a cork attached to a mounting plate. If a freezer that was originally built with a liner with a cork insert needs to be replaced, you will need to use the cork and mounting plate assembly as a replacement. The field service replacement PN is 141125215. This assembly includes the cork mounted on the mounting plate as well as the gasket and rivets required to install them. Units that only need a gasket replacement can retain the existing liner and replace only the gasket.

**Q:** Can I use Chart Connect 3000 with a TEC 3000 controller? What about using the TEC Connect program with a TEC 2000?

**A:** Both Chart Connect 3000 and TEC connect are compatible with TEC 2000 and TEC 3000's controllers. TEC Connect is the currently supported MVE controller interface and has increased functionality over Chart Connect 3000.

**Q:** The hot gas bypass temperature is not decreasing on a freezer that has just started filling, what should I do?

A: First verify that the supply source is connected properly and all supply valves are open. Verify that the gas bypass is enabled in the TEC 3000 menus, the bypass temperature set point is at a reasonable value (default value is -70 °C) and the gas bypass temperature probe is connected to the controller wiring harness (if the probe is not connected the bypass temperature will read "OPEN").

### ACCESSORIES

## Transfer Hose for Manual Discharge Device

As an alternative to the standard rigid spout on the manual discharge device, Chart offers a flexible, stainless steel transfer hose with either a small or large phase separator. The hose is 2 feet long with 1/4 inch diameter threads.



Transfer Hose for Discharge Device

Size	Part Number
Small	14044143
Large	14044151



Stainless Steel Transfer Hose in Use

8]glf]Vi hYX'Vm PrincetonCryo.com Cryogenic Delivery and Storage Systems

www.PrincetonCryo.com | Sales@PrincetonCryo.com | 800.232.2796