



MVE

Tech Tips



A monthly publication for the MVE Biological Products Distributors

MVE Replacement Kits Now Available

It has been requested by several of our distributors that kit assemblies for brass fittings used on XLC freezers are made available for purchase. These brass parts are vulnerable to cracking if plumbing assemblies are subjected to extreme vibrations, specifically the parts that connect to stainless. The kit numbers are as follows:

11081660	KIT XLC CABINET BRS PARTS
11081678	KIT 511/1211/1841 BRS. CONN.
11081686	KIT 810/1370/HE'S BRS. CONN.
11081694	KIT 230 GEN 2 BRS. CONN.

These kits make up the brass elbows and reducers used on the plumbing of XLC's and would be very beneficial to keep on hand when performing service calls.

Last January, MVE introduced the Cryo-claw. This handy device is used for retrieving and extracting vials, canes, straws or boxes that are stored in nitrogen freezers. Over time, the suction cups on this retrieving device may become brittle and break. Replacement cups are now available from MVE. To order a set you can contact Customer Service and request part number 11205221.

MVE No Longer Services Gordinier Electronics Automatic Controllers

Prior to 1993, MVE freezers used automatic controllers manufactured by Gordinier Electronics. These can be found on the A-4500, XLC-1110, and the XLC-440. MVE can no longer service these controllers or supply electronic components. Only the sensor assembly (part number 4613349) will still be available from MVE. This part can only be used if the freezer has full automatic capabilities.

If problems develop with the electronics on these systems and replacements parts or repair is required you can contact Gordinier Electronics at 810 778-0426. A return/repair program can be set up with them directly.

8]ghf]Vi hX'Vm



PrincetonCryo.com
Cryogenic Delivery and Storage Systems