



Tech Tips



September 2010

IATA Cryovials

Many varieties of vials are available to our customers today. Chart offers plastic, externally threaded Cryovials as an IATA approved accessories. Each part number is 100 vials.

Externally Threaded Cryovials

Description	Part Number
CRYOVIAL IATA PACKAGE 1.2 ML	10778810
CRYOVIAL IATA PACKAGE 2ML	10778828
CRYOVIAL IATA PACKAGE 3ML	10778836
CRYOVIAL IATA PACKAGE 4ML	10777438



Manual Discharge Device Transfer Hoses

Several filling accessories are available for the manual discharge device. The stainless steel transfer hoses are two feet long, have threads 1/4 inch in diameter, and are equipped with a choice of two phase separators.

Stainless Steel Transfer Hose

Phase Separator Size	Part Number
Small	14044143
Medium	14044151



Tec 2000 to TEC 3000 Upgrade Wiring, PN: 13319512

When wiring your stand alone TEC 3000 upgrade kit for MDC TEC 2000 controllers, please be mindful of the battery backup polarity. Also, be sure you are working from the correct instruction revision. The document should read PN 13319601 Rev A.

Distributed by



PrincetonCryo.com
Cryogenic Delivery and Storage Systems

www.PrincetonCryo.com | Sales@PrincetonCryo.com | 800.232.2796